# euro gine

Technology for Women's Healthcare

Rev. 02/2023

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## INFORMATION FOR THE USER

Read this leaflet carefully with your doctor and make sure you are told everything you need to know about this method at any time. Ask any questions you might have. This information is not intended to replace the explanation given by the health provider.

### IUD DESIGN, SIZE AND COMPOSITION

The intrauterine contraceptive device (IUD) with copper consists of a polyethylene structure with copper wire wrapped around its vertical shaft. A polyethylene thread is knotted to the lower part of this shaft and is used to remove the device.

Depending on the shape of the woman's uterus, the doctor can choose from up to 2 different sizes of IUD (depending on the model). At the same time, the doctor may choose between two shapes (**Eurogine® T 380, Gold T®** or **Ancora 375**; depending on the model) and solid copper wire with a silver or gold core (depending on the model).

**Normal size:** 31 mm wide x 33 mm long (**Eurogine® T 380, Gold T®**) / 20 mm wide x 35 mm long (**Ancora 375**)

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**Composition:** polyethylene, barium sulphate (radiopaque to x-rays), titanium oxide.

### MECHANISM OF ACTION AND POSSIBLE EFFECTS ON THE MENSTRUAL CYCLE

The IUD is inserted into the uterus to prevent pregnancy by influencing sperm motility and inhibiting their capacity to fertilise an egg before reaching the uterus.

During the first few months after insertion of the IUD occasional bleeding or spotting may occur between periods, which will remit over time (3-6 months), as will any cramp, abdominal and/or lower back pains similar to those suffered by some women during menstruation; menstrual bleeding may be more intense or prolonged than normal.

The pregnancy rate with an IUD is very low and similar to oral methods (especially for high load models with 375/380 mm<sup>2</sup> of active copper surface area) and its efficacy is maintained as of the first day after insertion.

### BENEFITS OF THE COPPER IUD

- High long-term efficacy and low rate of complications
- Does not interfere with sexual relations
- Does not interfere with breastfeeding
- Fertility is regained immediately after removal
- Does not affect ovulation
- No period of rest is required
- Does not have the typical side effects of medicated devices
- Does not contain hormones
- According to the WHO, it is the best contraceptive method in terms of cost-effectiveness and safety
- Can be used as emergency contraception
- According to studies, it may be a protective factor against the risk of developing cervical or endometrial cancer

### CONTRAINDICATIONS AND SPECIAL PRECAUTIONS

- Known or suspected pregnancy
- Puerperal sepsis
- Incipient Infection/Inflammation/PID
- Untreated acute cervicitis or vaginitis or other lower genital tract infection
- Post-partum endometritis or immediately after septic or second-trimester abortion
- Undiagnosed abnormal uterine bleeding
- Untreated malignant gynaecological disease
- Uterine anomalies
- Allergy to copper
- Allergy or hypersensitivity to any indicated component of the IUD including copper or any other trace elements present in copper
- Wilson's disease
- Strong susceptibility to infections
- Previous inserted IUD and not removed
- Severe dysmenorrhoea, unexplained vaginal bleeding
- Blood clotting disorders
- Cardiovascular or heart disease
- Anaemia/thalassaemia
- Endometriosis
- STD
- AIDS, antiretroviral therapy

#### Special precautions:

Consider the suitability of the method in the event of anaemia, dysmenorrhoea, and use of anticoagulants, non-steroidal anti-inflammatory drugs or corticosteroids.

Tampons should not be used during menstruation. If used, take care when removing them so as not to pull on the IUD removal threads. The same is so with menstrual cups.

### DRUG INTERACTIONS AND OTHER TREATMENTS

Tell the doctor prescribing drugs that you use a copper IUD. Reports have been published indicating that the contraceptive effect of IUDs seems to be reduced in patients receiving chronic treatment with non-steroidal anti-inflammatory drugs (particularly acetylsalicylic acid) and corticosteroids; nevertheless, experience has shown that the contraceptive protection is not reduced during short-term treatment of dysmenorrhoea with non-steroidal anti-inflammatory drugs.

High-frequency current electrotherapy (diathermy, shortwave or microwave) should be performed with caution in the vicinity of the abdomen and sacrum and it is important to follow the instructions of the manufacturer of the equipment used.

It has been shown that the use of magnetic resonance imaging (MRI) systems of up to 3 Tesla are safe for copper IUDs.

When receiving radiotherapy, the patient must inform the doctor that she is using an IUD.

### INSERTION AND REMOVAL PROCEDURES

Insertion/removal may cause some degree of pain and bleeding, fainting due to a vasovagal reaction or seizures in the case of epilepsy.

**Insertion:** The doctor will use a vaginal speculum to view the uterine cervix and clean it with antiseptic solution. The IUD is then inserted into the uterus via a thin, flexible tube which you will be able to feel. Finally, the doctor will cut the removal threads. The IUD can be inserted at any stage of the menstrual cycle; nevertheless, it should preferably take place during the last few days of menstruation to rule out the possibility of pregnancy and, furthermore, the cervical canal is dilated, and this simplifies insertion.

**Removal:** The best time for removal or replacement of the IUD is during menstruation because the cervical canal is slightly dilated, and pregnancy is ruled out. If it is removed halfway through the cycle and the user has had sexual intercourse during the previous week, there is a risk of pregnancy unless the IUD is replaced immediately.

After inserting a speculum to view the cervix, forceps are used to gently pull on the IUD threads until the device is fully removed.

It may be replaced immediately after removal under normal conditions. Unless another IUD is inserted, you are not protected against unwanted pregnancy.

### PROCEDURE FOR CHECKING THE PRESENCE OF THE IUD

At least once a month, and very especially immediately after each period, make sure that the IUD removal threads are still in place. Your health provider must teach you the procedure. Gently insert a clean finger into the vagina until you feel the strings close to the cervix. Do not pull on the threads. If you feel any change in their length, do not feel them or you can feel any other part of the IUD different from the threads, the IUD may not be in its place or may have been totally expelled. If you cannot feel the threads, take additional contraceptive measures -you may not be protected against unwanted pregnancy- and contact your doctor.

### RISKS AND COMPLICATIONS

**Perforation of the uterus or cervix:** this may occur at the time of insertion or subsequently, although the risk is greater when inserted during breastfeeding and for up to 36 weeks after birth. Other factors could be medical malpractice or abnormal uterine anatomy. Delay in the diagnosis of a perforation may reduce the contraceptive effectiveness or the IUD could migrate out of the uterine cavity. Perforation may reduce the contraceptive effectiveness.

**Sepsis:** acute infections or sepsis have been reported, including Group A Streptococcal Sepsis (GAS), after IUDs insertion. In some cases, a severe pain occurred within hours after insertion, followed by sepsis within days. As death from GAS is more likely if treatment is delayed, it is important to be aware of this rare but serious infections.

**Pelvic infections, Pelvic Inflammatory Disease and Endometritis:** if pelvic inflammatory disease (PID) or other infections occur, this usually happens in the 20 days following insertion and is often associated with an STD. It can interfere with fertility and increase the risk of ectopic pregnancy, to require hysterectomy, to cause sepsis or death. PID and Endometritis may be asymptomatic but still result in serious complications.

**Partial or total expulsion:** it mainly happens in the first year after insertion, especially in the first 3 months, and often with menstruation. The symptoms that suggest total or partial expulsion can include pain, bleeding and/or noticing that the threads are absent or longer than expected. Some expulsions are asymptomatic. Partial expulsion can reduce the device's contraceptive efficacy. A total expulsion prevents the contraceptive efficacy.

**Pregnancy:** the IUD is a very effective method of contraception, but it is not 100% effective. A pregnancy in the presence of an IUD is most likely to occur during the first year of use. It is recommended that the IUD should be removed, regardless of whether or not you wish to continue with the pregnancy. Continuing with an IUD in situ could lead to complications, the most common complication being miscarriage (including septic abortion).

**Ectopic pregnancy:** using an IUD does not increase the risk of an ectopic pregnancy, which is even lower than if no contraceptive method is used, although a pregnancy with an IUD in situ is more likely to be ectopic. It may require surgery and may result in loss of fertility.

**Not visible threads:** this may suggest retraction of the threads into the cervical canal or the uterus, total or partial expulsion of the IUD or perforation.

**Embedment:** the IUD may partially penetrate into the myometrium and make the removal difficult, even make the IUD break. In some cases, a surgical removal may be necessary.

**Breakage:** although very unlikely, the IUD could break during removal or in the uterine cavity. Its contraceptive efficacy could be compromised in some cases. Additional treatment may be needed.

**Allergy:** IUDs that include copper may cause allergic reactions.

### CONTACT YOUR DOCTOR IF:

- You miss a period or think that you could be pregnant.

- (en)
- During the first few days after insertion you may have slight bleeding, spotting and/or cramp and these symptoms do not abate with time or become worse.
  - You suffer persistent abdominal pain accompanied by unexplained fever or unusual and/or foul-smelling vaginal discharge, as you could have an infection.
  - Menstruation is painful and abnormal.
  - Your partner or you feel discomfort during sex. The IUD could have moved or caused cervical perforation, or the strings may be too long.
  - You cannot feel the removal threads during regular controls, as this means you could have expelled the IUD, or it could have moved or caused perforation. In these cases, the symptoms may be abnormal bleeding and/or abdominal pain.
  - You can feel the threads longer than usual or you can feel any other part of the IUD different from the threads during your regular controls. The IUD may be dislocated.
  - You notice the total or partial expulsion of the IUD.
  - You are prescribed long-term treatments with anticoagulants, non-steroidal anti-inflammatory drugs or corticosteroids.
  - You feel persistent itching or urticaria in the genital area, as this may be a possible allergy to copper or an infection of some type.
  - You may have been exposed to or you are diagnosed with a sexually transmitted disease (STD), HIV, AIDS or another gynaecological disease.
  - You have any doubts or concerns requiring clarification.

#### **OTHER WARNINGS**

- IUDs do not afford any protection against HIV infection or any other sexually transmitted disease (STD).
- Have regular check-ups according to your doctor's instructions.
- Ask for the traceability of the IUD and keep this information. It is included on the labels of the inner and outer packaging and must also be recorded by the doctor in your medical record.
- You must give your written consent to authorise insertion of the IUD.
- Eurogine can accept no liability if the instructions in this leaflet are not observed.

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